INSTRUCTION FOR USE
使用說明

A Guide to Soft Contact Lens Wear
軟式隱形眼鏡配戴指南

Table of Contents

English (EN):.....	2
繁體中文 (TC):	6

English (EN):

MiSight 1 day Soft Contact Lens

Instruction for Use

IMPORTANT: Please read carefully and keep this information for future use. This instruction for use is intended for the eye care practitioner, but should be made available to patients upon request. The eye care practitioner should provide the patient with the patient instructions that pertain to the patient's prescribed lens.

SYMBOLS KEY:

The following symbols may appear on the label or carton.

SYMBOL	DEFINITION
	Caution: this device to sale by or on the order of a licensed practitioner
	Caution / See Instructions for Wearers
	Use by Date (expiration date)
	Batch Code
	Sterile using Steam Heat
	Manufacturer
	Authorized representative in the European Community
	Do not use if package is damaged
	Consult instructions for use / consult electronic instructions for use
	Do not re-use
	Date of manufacture
	Medical device
	European Conformity Mark
L	Left
R	Right
SPH	Spherical power in diopters
BC	Base curve
DIA	Diameter
Made In	Made in statement / country of manufacture

DESCRIPTION

MiSight (omafilcon A) Contact Lenses are designed for myopia control and vision correction, indicated for single use daily disposable wear.

The lenses are made from a material containing 60% water and 40% omafilcon A, consisting of 2-hydroxy-ethylmethacrylate and 2-methacryloyloxyethyl phos-phorycholine polymers cross-linked with ethyleneglycol dimethacrylate. The lens material has a permanently fixed tint using Vat Blue 6, which is added to make the lens more

visible for handling.

Lenses parameters:

- o Diameter: 13.00 mm to 15.5 mm
- o Basic Curve: 8.00 mm to 9.50 mm
- o Center Thickness: 0.08 mm to 0.140 mm (varies with power)
- o Powers: -0.25D to -6.00D in 0.25D steps
-6.50D to -10.00D in 0.5D steps
- o Wearing Schedule Daily Wear
- o Replacement Schedule Daily Disposable
- o Lens Color Light Blue

The physical/optical properties of the lens are:

- o Refractive Index 1.395 ± 0.005
- o Light Transmittance > 90%
- o Water Content 60%
- o Surface Character hydrophilic
- o Oxygen Permeability 25 x10⁻¹¹ (cm²/sec)(ml O₂/ml x mmHg) (Polarographic FATT Method)

INDICATIONS FOR USE

MiSight (omafilcon A) Soft Contact Lenses may reduce the rate of myopia progression in children (6-18) and correct ametropia. Reduction of myopia progression was observed in children with a wearing time of 12 hours (8-16 hours) per day, 6.4 days (5-7) per week in a clinical study.

MiSight (omafilcon A) Soft (hydrophilic) Contact Lenses may provide improved comfort for contact lens wearers who experience mild discomfort or symptoms related to dryness during lens wear associated with Evaporative Tear Deficiency or from Aqueous Tear Deficiency (non-Sjogren's only).

ACTIONS

When placed on the cornea in its hydrated state, the **MiSight** (omafilcon A) Soft (Hydrophilic) Contact Lens acts as a refracting medium to focus light rays on the retina and to simultaneously provide an optical stimulus to slow the progression of myopia.

CONTRAINDICATIONS

Do not use the **MiSight** daily wear single use lens when any of the following conditions exist:

- o Acute and subacute inflammation or infection of the anterior chamber of the eye.
- o Any eye disease, injury, or abnormality that affects the cornea, conjunctiva, or eyelids.
- o Severe insufficiency of lacrimal secretion (dry eyes).
- o Corneal hypoesthesia (reduced corneal sensitivity), if not aphakic.
- o Any systemic disease that may affect the eye or be exaggerated by wearing contact lenses.
- o Allergic reactions of ocular surfaces or adnexa that may be induced or exaggerated by wearing contact lenses or use of contact lens solutions.
- o Any active corneal infection (bacterial, fungal, or viral).
- o If eyes become red or irritated.
- o The patient is unable to follow lens care regimen or unable to obtain assistance to do so.

WARNINGS

Patients were not studied who exceed the conditions characterized by any of the following diagnostic parameters:

- o Rose Bengal staining > 12 on a scale of 18
- o Fluorescein staining >12 on a scale of 15

- Meibomian gland dysfunction >3 on a scale of 0-4

Patients should be advised of the following warnings pertaining to contact lens wear:

- **PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYE.**
- Patients should be cautioned that proper use and care of the contact lenses and lens care products, including lens cases, are essential for the safe use of these products.
- It is essential that patients follow their eye care practitioner's directions and all labeling instructions for proper use of lenses and lens care products.
- Patients should fill their lens case with fresh solution every time they store their lenses, and never re-use solution. Additionally, they should clean and rinse their lens case between uses as recommended by their eye care practitioner.
- If patients experience eye discomfort, excessive tearing, vision changes, redness of the eye or other problems, they should be instructed to immediately remove their lenses and promptly contact their Eye Care Practitioner. It is recommended that contact lens wearers see their Eye Care Practitioner routinely as directed.
- Do not expose the contact lenses to water while wearing them.
- Water can harbor microorganisms that can lead to severe infection, vision loss or blindness. If the lenses have been submerged when swimming in pools, lakes, or oceans, discard them and replace them with a new pair.
- Ask eye care practitioner (professional) for recommendations about wearing the lenses during any activity involving water.
- Eye problems, including corneal ulcers, can develop rapidly and lead to **loss of vision**. Patient should be instructed to immediately remove the lenses, and promptly contact the eye care practitioner.
- Daily wear lenses are not indicated for overnight wear, and patients should be instructed not to wear lenses while sleeping. Clinical studies have shown that risks of serious adverse reactions is increased when these lenses are worn overnight.
- Single Use, Daily Disposable lenses are not intended for cleaning or re-use, and on removal should be discarded and a fresh pair used each day.
- Studies have shown that contact lens wearers who are smokers have a higher incidence of adverse reactions than non-smokers.
- MiSight lenses provide an optical correction that simultaneously presents one image "in-focus" and a second image "out-of-focus". Under certain circumstances (such as low light levels), this optical design can cause the following visual symptoms for some patients:
 - Reduced image contrast;
 - A ghost image (double image, with a mild second image seen);
 - Halos around bright lights; and glare around lights.
- Not all patients function equally well with this type of correction. This type of correction can create a vision compromised that may cause difficulties with certain visually-demanding tasks. Patients should exercise extra care if performing potentially hazardous activities. The overall risk of ulcerative keratitis may be reduced by carefully following directions for lens care, including cleaning the lens case.

PRECAUTIONS

Special Precautions for Eye Care Practitioners

- When selecting an appropriate lens design and parameters,

the eye care practitioner should consider all characteristics of the lens that can affect lens performance and ocular health, including oxygen permeability, wettability, central and peripheral thickness, and optic zone diameter.

- The potential impact of these factors on the patient's ocular health should be carefully weighed against the patient's need for refractive correction; therefore, the continuing ocular health of the patient and lens performance on the eye should be carefully monitored by the prescribing eye care practitioner.
- Fluorescein, a yellow dye, should not be used while the lenses are on the eyes. The lenses absorb the dye and become discolored. Whenever fluorescein is used in the eyes, the eyes should be flushed with a sterile saline solution that is recommended for in-eye use.
- Before leaving the eye care practitioner's office, the patient should be able to promptly remove the lenses or should have someone else available who can remove the lenses for him or her.
- Eye care practitioners should instruct the patient to remove the lenses immediately if the eye becomes red or irritated.

Eye care practitioners should carefully instruct patients about the following care regimen and safety precautions:

- Lenses prescribed on a daily wear single use wearing schedule should always be discarded when removed at the end of the wearing day.
- Do not use saliva or anything other than the recommended solutions for lubricating or wetting lenses.
- If the lens sticks (stops moving) on the eye, follow the recommended directions on Care for a Sticking Lens. The lens should move freely on the eye for the continued health of the eye. If non-movement of the lens continues, the patient should be instructed to **immediately** consult his or her eye care practitioner.
- Always wash and rinse hands before handling lenses. Do not get cosmetics, lotions, soaps, creams, deodorant, or sprays in the eyes or on the lenses. It is best to put on lenses before putting on makeup. Water-based cosmetics are less likely to damage lenses than oil-based products.
- Do not touch the contact lenses with the finger or hands if the hands are not free of foreign materials, as lens damage may occur.
- Carefully follow the handling, insertion, removal, cleaning, and wearing instructions in the Patient Instructions for **MiSight** contact lenses and those prescribed by the eye care practitioner.
- Never wear lenses beyond the period recommended by the eye care practitioner.
- If aerosol products such as hairspray are used while wearing lenses, exercise caution and keep eyes closed until the spray has settled.
- Always handle lenses gently and avoid dropping them.
- Avoid all harmful or irritating vapors and fumes while wearing lenses.
- Ask the eye care practitioner about wearing the lenses during sporting activities.
- Lenses prescribed on a daily disposable wearing schedule should always be discarded when removed at the end of the wearing day.
- Inform the doctor (health care practitioner) about being a contact lens wearer.
- Never use tweezers or other tools to remove lenses from the lens container unless specifically indicated for that use. Pour the lens into your hand.
- Do not touch the lens with fingernails.
- Always contact the eye care practitioner before using any

- medicine in the eyes.
- Always inform the employer of being a contact lens wearer. Some jobs may require use of eye protection equipment or may require that the patient not wear contact lenses.
- As with any contact lens, follow-up visits are necessary to assure the continuing health of the patient's eyes. The patient should be instructed as to a recommended follow-up schedule.

ADVERSE REACTIONS

The patient should be informed that the following problems may occur when wearing contact lens:

- Eyes stinging, burning, or itching (irritation), or other eye pain.
- Comfort is less than when the lens was first placed on the eye.
- Feeling that something is in the eye such as a foreign body or a scratched area.
- Excessive watering (tearing) of the eyes.
- Unusual eye secretions.
- Redness of the eyes.
- Reduced sharpness of vision (poor visual acuity).
- Blurred vision, rainbows, or halos around objects.
- Sensitivity to light (photophobia).
- Dry eyes.

If the patient notices any of the above, he or she should be instructed to:

- **Immediately remove the lenses.**
- If the discomfort or the problem stops, then look closely at the lens. If the lens is in some way damaged, **do not** put the lens back on the eye. If the problem continues, **do not** put the lens back on your eye; **immediately remove the lenses and consult the eye care practitioner.**
- When any of the above problems occur, a serious condition such as infection, corneal ulcer, neovascularization, or iritis may be present. The patient should be instructed to **keep the lens off the eye and seek immediate** professional identification of the problem and prompt treatment to avoid serious eye damage.

FITTING

Conventional methods of fitting contact lenses apply to all **MiSight** contact lenses.

WEARING SCHEDULE

For best results, it is recommended that the patient wears the lens for a minimum of 10 hours per day for at least 6 days per week. Daily wear lenses are not indicated for overnight wear, and patients should be instructed not to wear lenses while sleeping.

All **MiSight** lenses should be discarded and replaced with a new lens on a daily basis.

CARE FOR A DRIED OUT (DEHYDRATED) LENS

If any **MiSight** lens is exposed to air while off the eye, it may become dry and brittle. In this event, simply dispose of the lens and replace with a fresh one.

CARE FOR A STICKING (NONMOVING) LENS

If the lens sticks (stops moving or cannot be removed), the patient should be instructed to apply 2 to 3 drops of the recommended lubricating or rewetting solution directly to the eye and wait until the lens begins to move freely on the eye before removing it. If non-movement of the lens continues more than 5 minutes, the patient should immediately consult the eye care practitioner.

EMERGENCIES

The patient should be informed that if chemicals of any kind (household products, gardening solutions, laboratory chemicals, etc.) are splashed into the eyes, the patient should: **FLUSH THE EYES IMMEDIATELY WITH TAP WATER AND IMMEDIATELY CONTACT THE EYE CARE PRACTITIONER OR VISIT A HOSPITAL EMERGENCY ROOM WITHOUT DELAY.**

HOW SUPPLIED

Each lens is supplied sterile in a blister containing buffered isotonic saline solution. The blister is labeled with the base curve, diameter, dioptric power, manufacturing lot number, and expiration date of the lens.

DO NOT USE IF THE MISIGHT LENS IS BROKEN OR THE SEAL HAS BEEN DAMAGED

REPORTING OF ADVERSE REACTIONS

All serious adverse experiences and adverse reactions observed in patients wearing any **MiSight** contact lens or experienced with the lenses should be reported to:

[Local Responsible Person]

CooperVision (HK) Ltd.

Unit Nos. 1805-1806, Level 18, 909 Cheung Sha Wan Road,
Kowloon, Hong Kong
Telephone: (852) 37180699
Fax: (852) 24261177

[Manufacturer]

Coopervision, Inc.

711 North Road, Scottsville, New York 14546, USA

[Manufacturing Site]

1. CooperVision, Inc.
711 North Road, Scottsville, New York 14546, USA
2. CooperVision Manufacturing Ltd.
South Point, Hamble, Southampton, S031 4RF, United Kingdom
3. CooperVision Manufacturing Puerto Rico LLC,
500 Road 584, Lot 7, Amuelas Industrial Park, Juana Diaz, 00795 Puerto Rico
4. CooperVision Manufacturing Costa Rica, SRL
Zonal Franca El Coyol, Building 53, Alajuela, 20101 Costa Rica
5. CooperVision CL KFT.
Gorcsev Ivan Street 7, Building C, 2360 Gyal, Hungary

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1. CooperVision, Inc.
711 North Road, Scottsville, New York 14546, USA
2. CooperVision Manufacturing Ltd.
South Point, Hamble, Southampton, S031 4RF, United Kingdom
3. CooperVision Manufacturing Puerto Rico LLC,
500 Road 584, Lot 7, Amuelas Industrial Park, Juana Diaz, 00795
Puerto Rico
4. CooperVision Manufacturing Costa Rica, SRL
Zonal Franca El Coyol, Building 53, Alajuela, 20101 Costa Rica
5. CooperVision CL KFT
Gorcsev Ivan Street 7, Building C, 2360 Gyal, Hungary