



CooperVision®
Live Brightly.™

一樣 · 其實不一樣 clariti® 1 day 全面進化

2014

2019



- 1. 全新的鏡片邊緣設計大幅提升舒適度
- 2. 投資超過2億美元發展先進環保的生產技術
- 3. 優化製造技術確保高質素及精準的鏡片生產
- 4. 品質控制達至Biofinity系列標準



3X
高透氧

clariti® 1 day 高透氧隱形眼鏡已經全面進化，是時候為你的顧客升級轉用矽水凝膠隱形眼鏡

clariti® 1 day
see life with clariti®

比起一般水凝膠隱形眼鏡，clariti® 1 day 為角膜提供3倍高透氧，讓你顧客的眼睛更淨白明亮，體驗健康

新發佈

clariti® 繼續進化 最多功能 · 最齊度數¹



讓你的顧客得到更多


立即為你的散光顧客升級 clariti® 1 day **散光** 鏡片，一款最多度數組合的每日即棄矽水凝膠散光隱形眼鏡¹



clariti® 1 day **散光** 隱形眼鏡是：



每日即棄矽水凝膠散光鏡片中**最多度數組合**的產品¹



度數組合增加接近50%



clariti® 1 day toric has **AROUND-THE-CLOCK AXES****

線位每10度一跳為你照顧更多散光顧客的需要**

更多度數組合*：

單光度數	散光	線位
+0.00D to -6.00D (0.25D steps)	-0.75, -1.25, -1.75	10°, 20°, 30°, 40°, 50° , 60°, 70°, 80°, 90°, 100°, 110°, 120°, 130°, 140°, 150° , 160°, 170°, 180°
	-2.25	10°, 20°, 70°, 80° , 90°, 100°, 110° , 160°, 170°, 180°
-6.50D to -9.00D (0.50D steps)	-0.75, -1.25, -1.75	10°, 20°, 60°, 70°, 80°, 90°, 100°, 110°, 120°, 160°, 170°, 180°
	-2.25	10°, 20°, 90°, 160°, 170°, 180°
+0.25D to +4.00D (0.25D steps)	-0.75, -1.25, -1.75	10° , 20°, 70°, 80° , 90°, 100° , 110°, 160°, 170° , 180°

紫色的數字代表最新添加的線位

clariti® 1 day 散光鏡片為你的顧客提供更健康[†]、更貼身的體驗

* Prescription includes sphere power, cylinder, axis and ADD combinations. **Around-the-clock axes in 10° from plano to -6.00DS in -0.75DC, 1.25DC, and -1.75DC. References: 1. CVI data on file 2019. Based on number of published power options available across all soft lenses from Alcon, Johnson and Johnson, and Bausch + Lomb, UK and US April 2019.[†] With higher oxygen delivery than hydrogel materials, SiHy materials minimize or eliminate hypoxia-related complications during daily wear.